## INSTRUCTIONS FOR USE (100-test and 50-test kits)

### INTENDED USE

*SICKLEDEX*® is a qualitative solubility test kit for testing for the presence of sickling hemoglobins in human blood or sickle cell control material.

### SUMMARY AND PRINCIPLES<sup>1,2</sup>

Sickle cell disease is an inherited condition characterized by the presence of Hemoglobin S (Hb-S). Hb-S exists in a homozygous state (S/S) known as Sickle Cell Anemia or in a heterozygous state (A/S) known as Sickle Cell Trait. Homozygous individuals (S/S) commonly exhibit symptoms of severe hemolytic anemia and/or vascular occlusions. Heterozygous individuals (A/S) are usually asymptomatic. Hb-S may be present with other hemoglobins, such as Hemoglobin A, C or D, or with thalassemia, a condition that interferes with the synthesis of normal hemoglobin.

Under conditions of low oxygen tension, the heterozygous (A/S) form can cause erythrocytes to form the characteristic sickle-shaped tactoids. The formation of these irreversibly sickled red blood cells causes the onset of the acute symptoms. Detection of both the homozygous and heterozygous condition is important so high-risk individuals can be identified and their symptoms reduced.

Deoxygenated Hb-S is insoluble in the presence of a concentrated phosphate buffer solution and forms a turbid suspension that can be easily visualized. Normal Hemoglobin A and other hemoglobins remain in solution under these conditions. These different qualitative outcomes allow for the detection of sickle cell disease and its traits.

*SICKLEDEX* uses Saponin to lyse the red blood cells. Sodium Hydrosulfite then reduces the released hemoglobin. Reduced Hb-S is insoluble in the concentrated phosphate buffer and forms a cloudy, turbid suspension. Other sickling hemoglobin subtypes may also give a positive result.

### REAGENTS

*SICKLEDEX* Solubility Buffer is a 2.3M potassium phosphate buffer solution. The preservative is 0.1% 2-chloroacetamide. *SICKLEDEX* Reagent Powder vials contain Saponin and Sodium Hydrosulfite. Warning: Sodium Hydrosulfite is a flammable solid and strong reducing agent. Refer to SDS.

### PRECAUTIONS

- For In Vitro Diagnostic Use.
- Patient samples screened with this kit should be treated as if they are capable of transmitting infection and should be handled with appropriate precautions. Avoid contact with skin and mucous membranes.
- Reagent powder: Hydrosulfite: Harmful!

Risk (R): May cause fire. Harmful if swallowed. Contact with acid liberates toxic gases.

Safety (S): Keep container tightly closed and dry. In case of contact with eyes, rinse with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap. In case of fire, use dry powder. Never use water. Refer to SDS.

SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

Disposal of mixed buffer and reagent powder: Recommend disposal of mixed buffer/powder solutions in biohazardous waste bins.

### STORAGE AND STABILITY

*SICKLEDEX* is stable throughout the expiration date when stored tightly capped at 2°C to 30°C. **After opening and reconstituting the buffer solution, the working solubility buffer is stable for 45 days when stored tightly capped at 2°C to 10°C.** Allowing the bottle of reconstituted buffer to warm to room temperature (18°C to 30°C) may reduce its open-vial stability. DO NOT FREEZE. A slight sediment may form during storage. This will not interfere with test results.

### INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. The Saponin component in *SICKLEDEX* Solubility Buffer will cause the reagent to take on a yellow tint. This tint is normal and is not an indication of product deterioration. If the recovered test results are not as expected:

- Review control product package insert and the reagent test kit instructions.
- Check expiration date of the product on the vial. Discard outdated products.
- Repeat the run to confirm erroneous results; if the same result is obtained, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

### KIT CONTENTS

Materials Provided:

- SICKLEDEX* Solubility Buffer
- SICKLEDEX* Solubility Reagent Powder
- Dispenser caps for Solubility Buffer

Materials Required but not Supplied:

- 12 x 75mm disposable glass or plastic test tubes
- Test tube rack (paperboard or rack with lines)
- 20µl precision pipette with disposable tips or 20µl micropipettes
- Positive and negative Sickle-Chex® hemoglobin quality control samples are recommended. Sickle-Chex controls are available from Streck.

### REAGENT PREPARATION

The working solubility buffer must be prepared before screening can be performed.

- Bring buffer and reagent powder to room temperature before mixing.
- Add the contents of one vial of *SICKLEDEX* Reagent Powder to one bottle of *SICKLEDEX* Solubility Buffer.
- Place a white dispenser cap on the bottle of working solubility buffer. Dissolve the reagent powder completely with vigorous agitation.
- Record the reconstitution date and the reconstituted solution expiry date in the space provided on the solubility buffer bottle.
- Store the working solubility buffer tightly capped at 2°C to 10°C when not in use. Reconstituted buffer must be used within 45 days.

### BLOOD SAMPLE COLLECTION

- Fresh blood samples may be collected from a finger puncture.
- Use anticoagulated whole blood, packed cells, blood bank segments containing whole blood or packed cells with additive solutions. Never use clotted blood.
- Streck *SICKLEDEX* Solubility kits have been approved for use with the following commercial blood collection anticoagulants: EDTA, ACD, CPD, CPDA-1, CPDA-2, CP2D, Sodium Citrate and Heparin. Additives such as AS-1 and AS-3 used for blood banking do not affect *SICKLEDEX* results.
- Blood samples stored at 1°C to 10°C for up to 45 days may be used for testing.

### INSTRUCTIONS FOR USE

- Dispense 2.0ml of working *SICKLEDEX* Solubility Buffer into a 12 x 75mm disposable glass or plastic test tube. Return working solubility buffer to 2°C to 10°C immediately after use. Allow the working solution in the test tubes to warm to room temperature (18°C to 30°C) for a minimum of 10 minutes before use. The use of reagents below room temperature can give false results.
- Add 20µl of whole blood or 10µl of packed red blood cells to the test tube. When running control samples add 20µl of control. If the whole blood hematocrit is ≤15%, centrifuge the sample for 5-10 minutes at 1200 rpm. Pipet 10µl of the packed cell volume from the bottom of the tube and add it to the *SICKLEDEX* Solubility Buffer test tube.
- Mix the contents of the test tube thoroughly by swirling the tube several times. Place the test tube in the test tube rack. Optionally, secure the top of the tube with laboratory film and invert the tube several times as an effective mixing technique.
- Allow the sample to stand at room temperature (18°C to 30°C) for at least 6 minutes. Observe the sample for turbidity. Read results between 6 and 60 minutes.

### DIRECTIONS FOR USE OF 20µl PLASTIC MICROPIPETTE (STRECK PART #217663)

Never squeeze the plastic micropipette tube while sampling. Filling is automatic.

- Hold the tube horizontally, and touch the tip of the tube to the sample. Capillary action will automatically draw the sample to the fill line and stop.
- Dispense blood in the micropipette within 30 seconds to prevent clotting of the sample.
- To expel the sample, align the tip of the tube with the sample target and squeeze the bulb. If a sample won't expel, verify that the tube has been filled all the way. If necessary, touch the tip to the sample again and allow it to fill completely. Then align the tip with the sample target and squeeze the bulb to deliver the sample.


### LIMITATIONS<sup>1,2</sup>

- False positives may occur in patients with erythrocytosis, hyperglobulinemia, extreme leukocytosis or hyperlipidemia. Coarse flocculation may occur in these samples due to elevated levels of total serum protein. These patient samples may be washed in normal physiologic saline, centrifuged and 10µl of the packed cells used for testing.
- False positives or false negatives may occur in patients with severe anemia.
- The limit of detection is 15% Hgb-S with a total hemoglobin of 12g/dL
- 100% Sensitivity and 100% Specificity verified by testing 166 samples (137 negative and 29 positive).
- False negatives may occur in infants under six months of age due to elevated levels of Hemoglobin F.
- False positives or false negatives may occur in patients with a recent blood transfusion.
- Positive results may occur in patients with some rare sickling hemoglobin subtypes such as Hemoglobin C Harlem or Hemoglobin C Georgetown.
- SICKLEDEX* is a qualitative screening procedure and does not differentiate between Sickle Cell Disease (S/S) and Sickle Cell Trait (A/S). All positive test results should be further evaluated by hemoglobin electrophoresis, when used for patient testing. This does not apply to blood donor screening tests.
- SICKLEDEX* is compatible with RBC units that have been glycerolized and frozen, and subsequently thawed and deglycerolized, for the purpose of extending the storage life of an RBC unit. Each laboratory should validate their approved glycerolization protocol with *SICKLEDEX* prior to implementation.
- It should be noted that test time reads for some non-patient samples (proficiency or control material) may not be evident at 6 minutes and may take up to 60 minutes to be resolved.

### EXPECTED RESULTS

- The reaction is read macroscopically by looking through the test tubes at black lines of the test tube rack.
- A POSITIVE test for sickling hemoglobin is indicated by a cloudy, turbid suspension through which the black tube rack lines are NOT VISIBLE.
- A NEGATIVE test for sickling hemoglobin is indicated by a transparent suspension through which the black tube rack lines are CLEARLY VISIBLE.

Note: Sickle cell controls are designed to verify the activity of the reagent. Negative patient results may not clear as quickly or as completely as the control.

### REFERENCES

- McPherson R.A., Pincus M.R. Henry's clinical diagnosis and management by laboratory methods - 21st Edition. 2007.
- Miale J.B. Laboratory medicine hematology - Fourth Edition. 1972.

### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

### Rx Only

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

Streck  
7002 S. 109 Street, La Vista, NE 68128 USA

MediMask® Europe  
11, rue Emile Zola, BP 2332  
38033 Grenoble Cedex 2, France

350430-27  
2024-05

Instructions for use can be obtained at [streck.com](http://streck.com) in the product's Resources section, by calling 402-333-1982, or by calling your local supplier.

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