INSTRUCTIONS FOR USE

INTENDED USE

STaK-Chex® Plus Retics is an assayed whole blood control for evaluating the accuracy and precision of automated procedures that measure blood cell parameters.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation protocols.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument or manual methods, this whole blood control will provide values within the expected range indicated on the assay sheet. STaK-Chex Plus Retics is for use by laboratory professionals.

REAGENTS

These whole blood reagents may contain any or all of the following: stabilized human red blood cells, a nucleated red blood cell analog, a white blood cell component consisting of human or non-human analogs and a platelet component consisting of a human or non-human analog in a preservative medium. No additional materials are provided or required.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed of in general waste, but should be disposed of with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

Unopened STaK-Chex Plus Retics is stable 105 days when stored at 2 °C to 10 °C. After initial opening, STaK-Chex Plus Retics is stable 14 days when at 2 °C to 10 °C.

INDICATION OF PRODUCT DETERIORATION


Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. STaK-Chex Plus Retics should only be used with the instruments listed on the assay.
2. Remove vials of STaK-Chex Plus Retics from the refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. To mix: **(Do NOT mix mechanically or vortex.)**
For a video demonstration, visit streck.com/mixing.

- a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
 - d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
4. Refer to the instrument manual for analyzing control material.
 5. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and returning to refrigeration. Report any serious incident to Streck and appropriate regulatory entities including the competent member state in which the user and/or patient is established as applicable.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. If the product is value assigned for MCV and hematocrit, assay values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are

not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2, it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

'R' alarms may occur with control materials on automated systems. This will not affect the validity of the results. The alarms may be disregarded if the control is performing within the assay ranges.

Trending in the MCV parameter over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

Scanning the assay barcode will fill parameters in default instrument order. Confirm assigned target values are applied to correct parameters; if not, adjust accordingly.

CONTROL FILE SETUP PROCEDURES FOR BECKMAN COULTER® UniCel® DxH™ 600/800/900/690T

UniCel DxH 600/800/900/690T – QC SETUP

Combined CBC/Diff/Retic files

1. Select: **Menu/Setup/Quality Control**

2. Select: **New Patient Control**

3. Confirm the cursor selection is in the Specimen ID field: **Scan CBC/DIFF/Retics barcode on vial**

4. Choose Type: **CDR**

5. Choose Level

6. Select: **Create Control**

7. Confirm the cursor selection is in the Control ID field: **Scan 2-D barcode**

8. Select: **Save**

9. Repeat steps 2-8 for additional Control Levels

UniCel DxH 600/800/900/690T – RUNNING CONTROLS

Combined CBC/Diff/Retic

1. Place STKR vials in cassette with CBC/Diff/Retic barcode facing out.

2. Place cassette into output buffer

3. Data will automatically transfer to appropriate "Patient Control" file.

REFERENCES

1. Clinical Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers. Approved Standard – Second Edition.

QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional Information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

All product names, logos, brands, and marks are property of their respective owners.

See streck.com/patents for patents that may be applicable to this product.

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