## INSTRUCTIONS FOR USE

### INTENDED USE

CD-Chex Select® is intended to be used as a quality control material for evaluating monoclonal antibody binding by flow cytometry. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry they provide a reference for normal and abnormal peripheral blood leukocytes. CD-Chex Select is designed for use on BD® Biosciences and Beckman Coulter® flow cytometry systems. **This product and the markers provided on the assay have not been cleared by the U.S. Food and Drug Administration for In Vitro Diagnostic use. This product and the values provided are For Research Use Only. Not for use in diagnostic procedures.**

### SUMMARY/PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes present in a whole blood sample. Quantitation of peripheral blood leukocyte subsets has been shown to be useful in research and clinical applications including the evaluation of autoimmune and oncologic disorders. CD-Chex Select is a positive procedure assayed control used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance. CD-Chex Select mimics whole blood samples by possessing surface antigens and intracellular antigens detectable with fluorescent monoclonal antibodies by flow cytometry.

Subsets of leukocytes can be distinguished on the basis of light scatter properties as well as the presence of cell surface antigens. CD-Chex Select is designed to represent normal and abnormal peripheral blood leukocytes where lymphocytes, monocytes, and granulocytes are distinguishable on the basis of light scatter properties.

### REAGENTS

CD-Chex Select contains stabilized human blood in a preservative medium.

### PRECAUTIONS

1. CD-Chex Select is For Research Use Only. Not for use in diagnostic procedures.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates the use of the product.
5. CD-Chex Select should not be used as a calibrator.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

CD-Chex Select is stable through the expiration date when stored at 2 °C to 10 °C. After initial opening, CD-Chex Select is stable for 30 days when stored at 2 °C to 10 °C. DO NOT FREEZE.


### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of CD-Chex Select. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.
4. Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.

### INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):  
For a video demonstration, visit [streck.com/mixing](http://streck.com/mixing).  
a. Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- b. Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.


- c. Aliquot immediately after mixing.
- d. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.

4. Return control reagent to refrigeration immediately after sampling to ensure maximum open-vial stability.
5. Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.
6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
8. Analyze by flow cytometry using your laboratory's established protocol.

### GATING

The most common gating strategy in white blood cell assessment is to gate on the leukocyte subsets and then determine the CD marker percent positivity using a negative staining control<sup>1,2</sup>. Leukocytes are generally located on a FSC/SSC plot or a CD45/SSC plot; although other gating strategies can be employed<sup>1,2</sup>.

**FIGURE 1. Gating methods used to obtain reference values for CD-Chex Select.**

In flow cytometry, the lymphocyte, monocyte, and granulocyte populations can routinely be identified using a FSC/SSC plot (A and B). For CD-Chex Select, CD marker percent positivity for each leukocyte subset is determined using a negative staining control (C to N).

### INTRACELLULAR MARKER INSTRUCTIONS FOR USE

Follow steps 1-3 above. CD-Chex Select is a stabilized blood material. Therefore, the fixation step (Reagent 1 or Reagent A) used prior to permeabilization in commercially available intracellular staining kits is not required. Use of the fixation reagent will result in sub-optimal recovery.

### LIMITATIONS

1. Stabilized material is considered nonviable and not compatible with viability dyes and kits.

## EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers<sup>1,2</sup>. The assay values are obtained using common flow cytometry reagents. See the assay for limitations or specific instructions for reagents.

The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

## REFERENCES

1. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematolymphoid cells. Approved Guideline - Second Edition.
2. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.

## QUALITY CONTROL

Streck offers *STATS*<sup>®</sup>, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at [streck.com](https://streck.com).

## ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [streck.com](https://streck.com).

## GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [streck.com](https://streck.com).

All product names, logos, brands, and marks are property of their respective owners.

See [streck.com/patents](https://streck.com/patents) for patents that may be applicable to this product.

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