## INSTRUCTIONS FOR USE

### INTENDED USE

CD-Chex TdT Plus® is intended to be used as a quality control material for evaluating intracellular and surface antigens, including TdT (Terminal deoxynucleotidyl transferase), CD1a, CD34 and cytoplasmic CD3, with monoclonal antibody binding by flow cytometry. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry, they provide a reference value for abnormal cells found in certain types of hematopoietic neoplasms. CD-Chex TdT Plus is designed for use on BD® Biosciences and Beckman Coulter® flow cytometry systems. **This product and the markers provided on the assay have not been cleared by the U.S. Food and Drug Administration for In Vitro Diagnostic use. This product and the values provided are For Research Use Only. Not for use in diagnostic procedures.**

### SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry is a complex process. CD-Chex TdT Plus is designed to represent abnormal peripheral blood leukocytes consistent with a hematolymphoid neoplastic patient sample<sup>1,2,3</sup>. CD-Chex TdT Plus mimics whole blood samples by possessing surface antigens and intracellular antigens detectable with fluorescent monoclonal antibodies by flow cytometry. Stable antigens on CD-Chex TdT Plus include surface CD1a, CD2, CD4, CD5, CD7, CD8, CD34, intracellular CD3 and nuclear TdT. CD-Chex TdT Plus is a positive procedural assayed control used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance to provide consistent and reliable quality control measurements<sup>4</sup>.

### REAGENTS

CD-Chex TdT Plus contains stabilized human blood and cells of human origin in a preservative medium.

### PRECAUTIONS

1. CD-Chex TdT Plus is For Research Use Only. Not for use in diagnostic procedures.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. CD-Chex TdT Plus should not be used as calibrator.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

### STORAGE AND STABILITY

CD-Chex TdT Plus is stable through the expiration date when stored at 2 °C to 10 °C. After initial opening, CD-Chex TdT Plus is stable for 30 days when stored at 2 °C to 10 °C. DO NOT FREEZE.

### INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:


1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the CD-Chex TdT Plus. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.
4. Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.

### INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):

**For a video demonstration, visit [streck.com/mixing](https://streck.com/mixing).**

- a. Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- b. Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.

- c. Aliquot immediately after mixing.
- d. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.


4. Return control reagent to refrigeration immediately after sampling to ensure maximum open-vial stability.
5. Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.
 

Note: A negative staining control is recommended due to heterogeneous/dim expression of selected assayed parameters.
6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
8. Analyze by flow cytometry following the gating strategy outlined in the Gating section.

### GATING

The most common gating strategy used in neoplastic cell assessment is to gate on the abnormal cells and then determine the CD marker percent positivity using a negative staining control<sup>1,4</sup>. Abnormal cells are

generally located on a FSC/SSC plot or a CD45/SSC plot; although other gating strategies can be employed<sup>1-4</sup>. Debris should be excluded from the cell gate to obtain percent recovery values within assay range.

**FIGURE 1. Gating methods used to identify abnormal cells in CD-Chex TdT Plus.**

In flow cytometry, abnormal cells can routinely be identified by using a FSC/SSC (A and B) or a CD45/SSC (C and D). Optimal gates will remove debris for evaluation of percent recovery. Gates for percent recovery of each marker (E and F) were determined from unstained samples.

### INTRACELLULAR MARKER INSTRUCTIONS FOR USE

Follow steps 1-3 above. CD-Chex TdT Plus is stabilized cellular material. Therefore, the fixation step (Reagent 1 or Reagent A) used prior to permeabilization in commercially available intracellular staining kits is not required. Use of the fixation reagent will result in sub-optimal recovery.

### LIMITATIONS

1. CD-Chex TdT Plus is not intended to be used in ISHAGE protocols for CD34 enumeration. Staining characteristics are not equivalent to the phenotypic properties of human hematopoietic progenitor cells<sup>5,6</sup>.
2. Stabilized material is considered nonviable and not compatible with viability dyes and kits.

### EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers<sup>4</sup>. The assay values are obtained using common flow cytometry reagents. See the assay for additional limitations or specific instructions for reagents.

Use of alternative gating strategies not specified in these instructions may result in values outside the published range. The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

### REFERENCES

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematolymphoid cells. Approved Guideline - Second Edition.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

### QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

### ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

### GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

All product names, logos, brands, and marks are property of their respective owners.

See [streck.com/patents](https://streck.com/patents) for patents that may be applicable to this product.

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### RÉFÉRENCES

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematomalymphoid cells. Approved Guideline - Second Edition.

5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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## QUELLENANGABEN

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hemolymphoid cells. Approved Guideline - Second Edition.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.

- Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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## BIBLIOGRAFIA

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hemolympoid cells. Approved Guideline - Second Edition.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematology cells. Approved Guideline - Second Edition.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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## REFERENSER

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hemolymphoid cells. Approved Guideline - Second Edition.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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