INSTRUCTIONS FOR USE

INTENDED USE

Para 12[®] Plus is an assayed hematology control for evaluating the accuracy and precision of hematology instruments that provide a white blood cell differential.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation or manual protocols.

Use of stabilized cell preparations for controlling hematology instrumentation is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, Para 12 Plus will provide values within the expected range indicated on the assay sheet. Para 12 Plus is for use by laboratory professionals.

REAGENTS

This whole blood reagent may contain any or all of the following: stabilized human or mammalian red blood cells, human, mammalian or simulated white blood cells and a platelet component in a preservative medium. No additional materials are provided or required.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

Para 12 Plus is stable through the expiration date when stored unopened at 2 °C to 10 °C. After initial opening, Para 12 Plus is stable for 7 days when stored at 2 °C to 10 °C.

INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of Para 12 Plus. Discard outdated products.
3. Assay an unopened vial of Para 12 Plus. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.


INSTRUCTIONS FOR ONLINE ASSAY FILE DOWNLOAD

1. Insert a blank USB flash drive or 3½-inch floppy disk into the computer.
2. On the appropriate product page at streck.com, click the link for the assay disk file with the corresponding lot number you are using.
3. When prompted, Click "Open" to download and open the file.
4. Once the file is opened, click "File" → "Extract All" or "Extract All..."
Note: Windows[®] 7 Users - Be sure the box reading "Show extracted files when complete" is unchecked.
5. Follow the "Extraction Wizard" (if applicable). Click on the "Browse" button to select an extract destination. Select the flash drive or floppy that was inserted in Step 1.
6. After the files have been extracted, safely remove the flash drive or floppy disk from the computer and insert into the appropriate instrument.
7. Refer to instrument operator's manual for assay file set up.

INSTRUCTIONS FOR USE

Important Notes:

- ▶ The vial label barcode is only for use on the CELL-DYN Ruby[®].
 - ▶ CELL-DYN Enzymatic Cleaner Concentrate solution must be thoroughly removed from the analyzer Wash Block and/or Probe before running control products.
 - ▶ Para 12 Plus should only be used with the instruments listed on the assay.
1. Remove vials of control from the refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
 2. To mix: (**Do NOT mix mechanically or vortex.**)
For a video demonstration, visit streck.com/mixing.
 - a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.

d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).

3. Refer to the instrument manual for the system in use for analyzing control materials.
4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both vial and cap before replacing cap and returning to refrigeration. Report any serious incident to Streck and appropriate regulatory entities including the competent authority of the member state in which the user and/or patient is established as applicable.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with Para 12 Plus. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS AND THEIR DERIVATION

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. MCV and hematocrit values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2¹ it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

REFERENCES

1. Clinical Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers; Approved Standard – Second Edition.

QUALITY CONTROL PROGRAM

Streck offers STATS[®], an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

CHANGES FROM PREVIOUS VERSION

Updated to meet requirements of Regulation (EU) 2017/746.

See streck.com/patents for patents that may be applicable to this product.

All product names, logos, brands, and marks are property of their respective owners.

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